interscience

BagMixer® 400 W

WINDOW DOOR LIFETIME WARRANTY | SHOCK ABSORBERS LIFETIME WARRANTY | 3 YEAR WARRANTY

400 mL lab blender

Ref. 022 000 (100-127 V~ / 220-240 V~)
Ref. 022 230 (220-240 V~)



Standards:

ISO 7218 | ISO 6887-1 | FDA BAM Bacteriological Analytical Manual

Manufactured under:

RoHS 2002/95/EC | WEEE 2002/96/EC | CE | UK CA | MADE IN FRANCE

BagMixer® 400 W blends efficiently any kind of sample in any kind of bag, without any risk of contamination. The glass door allows a quick control of the blending in progress. BagMixer® 400 W is equipped with Click & Clean® (patented), a system of detachable paddles that facilitates cleaning by accessing the blending chamber entirely.

TECHNICAL SPECIFICATIONS

- All stainless steel 304L
- Window door
- Q-Tight® closing: tight closing for no more leaks from the top of the bag
- Fully opening door (270°)
- Security drip tray

- Useful volume: 50 - 400 mL
- Fixed speed (8 strokes/s)
- Variable blending time (30 - 210 s or ∞)

- Voltage and frequency
  Ref. 022 000 : 110-127 V~ / 220-240 V~ 50-60 Hz
  Ref. 022 230 : 220-240 V~ 50-60 Hz
- Fuses
  Ref. 022 000 : T 3.15 A L 250 V / T1.25A L 250 V 5x20 mm
  Ref. 022 230 : T 1.25 A L 250 V 5x20 mm
- Max power: 135 W

- Protection index: 21
- Dimensions (w x d x h): 26.8 x 44.1 x 28.3 cm, weight: 16.6 kg
- Box dimensions (w x d x h): 49 x 34.5 x 37 cm, gross weight: 19.7 kg
- 3-year warranty
- Spare parts availability: 10 years
- Shock absorbers lifetime warranty

BagMixer® 400 W_FT_06/22 Pictures and information are not contractually binding. BagMixer® is trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

interscience in the world

PARIS
Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT
Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON
Phone: +1 781 937 0007
sales.usa@interscience.com

SHANGHAI
电话: +86 178 2123 6642
sales.china@interscience.cn

SINGAPORE
Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO
Phone: +81 3 6712 9715
sales.japan@interscience.com

- Window door lifetime warranty
- Availability of spare parts: 10 years

- In compliance with ISO 7218, ISO 6887-1 and FDA-BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE
- Designed and made in France

DELIVERED WITH

- 1 security drip tray
- 1 power cable
- 1 user's manual
- 1 certificate of conformity / warranty card
- 1 free pack of blender bags

BagMixer® 400 W_FT_06/22 Pictures and information are not contractually binding. BagMixer® is trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

interscience in the world

PARIS
Phone: +33 1 34 62 62 61
info@interscience.com

FRANKFURT
Phone: +49 611 7238 7770
sales.germany@interscience.com

BOSTON
Phone: +1 781 937 0007
sales.usa@interscience.com

SHANGHAI
电话: +86 178 2123 6642
sales.china@interscience.cn

SINGAPORE
Phone: +65 6977 7232
sales.asia@interscience.com

TOKYO
Phone: +81 3 6712 9715
sales.japan@interscience.com