interscience

BagPipet®

Fixed-volume pipet & straw tips

Ref 251 091 (BagPipet®)



Standards:

ISO 7218 | ISO | ISO 4833-2 | FDA BAM Bacteriological Analytical Manual

Manufactured under:

MADE IN FRANCE

BagPipet® is an autoclavable pipet with fixed volume (0.1 mL / 0.9 mL / 1 mL). Used with BagTips® a single use rigid straw tip, available in 2 sizes (Regular 19 cm / Jumbo 24 cm), BagPipet® intakes any product, even very viscous products: milk, cream, yoghourt, cosmetics...The long BagTips® reaches the bottom of any bag.

TECHNICAL SPECIFICATIONS

- Fixed volume pipet
- Ejection of disposable straw
- Light and ergonomic
- Autoclavable
- Easy calibration
- 3 settings: 0.1 mL / 0.9 mL / 1 mL
- Calibrated according to the GLP “Good Laboratory Practice”
- Also available: volume push button
    - 100 µL (Ref 251 001)
    - 900 µL (Ref 251 009)
    - 1000 µL (Ref 251 010)
    - 100 à 1000 µL (Ref 251 080) Variable
- Dimensions : 6 x 18,5 x 3 cm, net weight: 0,149 kg
- Box (w x d x h): 30 x 10 x 6 cm, gross weight: 250 g
- 1 year warranty
- Availability of spare parts: 10 years
- In compliance with: ISO 7218, ISO 6887, ISO 4833-2 and FDA BAM (Bacteriological Analytical Manual)
- Designed and made in France

DELIVERED WITH

- 1 user’s manual
- 1 pouch BagTips® Regular
- 1 holder for key
- 1 volume button 1000 µL (mounted on the pipet)
- 2 volume buttons (900 µL, 100 µL)
- 2 O-Rings
- 1 Allen wrench
- Oil

ACCESSORIES

- Ref. 252 019 BagTips® Regular : Straw tip (5 to 12 kGy gamma ray treated) for BagPipet.
  Nominal volume: 1.1 mL - Length: 190 mm – Quantity per box: 1000
  Box (w x d x h): 27,5 x 23 x 11 cm, gross weight: 1,2 kg
- Ref. 252 024 BagTips® Jumbo : Straw tip (5 to 12 kGy gamma ray treated) for BagPipet.
  Nominal volume: 2,2 mL – Length: 240 mm - Quantity per box: 1000
  Box (w x d x h): 27,5 x 23 x 11 cm, gross weight: 1,5 kg

BagPipet® BagTips®_FT_0920 Pictures and information are not contractually binding. BagPipet® and BagTips® is trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

interscience in the world

PARIS Phone: +33 1 34 62 62 61 info@interscience.com
FRANKFURT Phone: +49 611 7238 7770 sales.germany@interscience.com
BOSTON Phone: +1 781 937 0007 sales.usa@interscience.com
SHANGHAI 电话: +86 178 2123 6642 sales.china@interscience.cn
SINGAPORE Phone: +65 6977 7232 sales.asia@interscience.com
TOKYO Phone: +81 3 6712 9715 sales.japan@interscience.com